

# NEWSLETTER 5/2009

FOR MEASURABLE SUCCESS

**SOEHNLE**  
PROFESSIONAL

Are you the wrong contact person for that Newsletter?  
Or would you like to have someone in your company to receive the Newsletter, too?  
Please email [us](#) the name and email of the contact person.

## What's new?


- \* **New product range of laboratory scales for precise measurements**  
**Available September 2009!**

Soehnle Professional introduces a new product range of laboratory scales and devices. We now offer precision and analytical scales. With their different functionalities and a modern design, they supplement Soehnle Professional's product range.

All scales are **approved either Class I or II**. A selection between three function levels - **Standard, Basic and Pro** - allows for finding the right product. Independent of the version, all scales offer:

- \* Backlighted display
- \* Functions: Piece counting, percentage weighing, checking, and underfloor weighing
- \* Internal calibration
- \* Serial interface RS232 as standard
- \* Connection of second indicator possible

### Precision scales, approved Class II



- \* **New product:**  
**Analogue/Digital Converter 3005**  
**Available September 2009!**

The **A/D Converter Professional 3005** enables the experienced user and scale manufacturer the installation of strain gauge load cell systems, which do not require an indicator.

The measuring data is provided for electronic processing via a selection of optional interfaces such as RS232, USB, Ethernet, Bluetooth™ and more. This allows you to easily connect scales to data processing systems. Up to six load cells in one scale are possible to connect.

Two versions are available:

The Basic version (3005.01.001) has **one measuring point for strain gauge load cells and one serial interface RS232 as standard**. The Profibus version (3005.05.001) offers additionally a Profibus.

Available functions are weighing, taring, switching, and dosing. With an internal division of 16 million digits and a measuring value rate of 80 per second, it is a perfect match for various industries such as plant construction, process technology, or even in scientific applications.

The system is easily configured via a PC software program, available in German, English and French.

If you are interested, please feel free to contact [us](#).

## Analytical scales, approved Class I




Our new range provides moisture analysers that are measuring devices especially designed for determining the moisture content of relatively small samples of different materials. The sample materials are first weighed and then dried. The difference of the initial weight and those after the drying process determine the moisture. The application area of moisture analysers: from the use in laboratories to industrial and commercial applications.

## Moisture analysers



We will inform you about our whole range of laboratory scales in the next weeks including a brochure with all details.



## Curios?

### MEDICA 2009

Soehnle Professional will exhibit again on the world's largest medical fair - **MEDICA 2009!**



18<sup>th</sup> - 21<sup>st</sup>  
November 2009

Meet us in Dusseldorf, Germany in hall 12 / stand E31 to learn about our latest products! We look forward meeting you!

Dusseldorf,  
Germany

[www.medica.com](http://www.medica.com)

Valid until 30<sup>th</sup> September 2009!

Please keep in mind our **special offer** in June/July 2009 for the **Compact scale 7755.12/15** and the **Baby Scale 8320!**



The Special offer is also available with more details [here](#).



## Contact:

You missed a Newsletter?

Follow the link <http://www.soehnle-professional.com/en/company/news/> and find the archive of our Newsletters.

If you have any questions or require further information, please do not hesitate to contact [us!](#)

Soehnle Professional GmbH & Co. KG  
Manfred-von-Ardenne-Allee 12  
71522 Backnang  
Germany  
Tel.: +49 - 7191 - 3453 - 272  
Fax: +49 - 7191 - 3453 - 211  
Email: [export@soehnle-professional.com](mailto:export@soehnle-professional.com)

[www.soehnle-professional.com](http://www.soehnle-professional.com)