

PRECISION SCALE

PRODUCT NUMBER 9231

Trade, commerce and industry



- > Weighing range: 120g, division 0.001g
- > Weighing range: 1,200g, division 0.01g
- > Back-lighted LCD display
- > Integrated accu
- > Calibration weight included

PRECISION SCALE

PRODUCT NUMBER 9231

Features



Weighing, counting, taring and percentage weighing. Precision scale with a very high resolution of 120,000 divisions. A highly precise measuring instrument for laboratories, recipes, checking and piece counting. The enclosed calibration weight guarantees accuracy at all times.

Technical data

- > Weighing platform made of stainless steel (ø 120 mm)
- > Backlighted LCD display with 13mm high digits
- > Bubble level
- > Interface RS232
- > Calibration weight included in delivery scope
- > Dust and splash proof IP54
- > Power supply 230V, integrated accu as standard
- > Battery saving mode

Functions

- > Weighing and taring
- > Counting
- > Percentage weighing
- > Plus-minus weighing
- > Unit switchover kg/lb

Optional

- > Wind shield (order number 2500.03.002) made of Plexiglas, sides pluggable as single parts, dimensions 182 x 145 x 100mm
- > Data cable (order number 2550.03.014)
- > Printer data cable (order number 2550.03.015)



Wind shield

Ordering information

Order number	Weighing range	Division	Linearity	Dimensions	Weight
	g	g	g	mm	kg
9231.01.001	120	0.001	± 0.01	180 x 230 x 65	1.2
9231.02.001	1,200	0.01	± 0.1	180 x 230 x 65	1.2

Your Soehnle Professional dealer:

www.soehnle-professional.com