

ADIPOSIITY SCALE

PRODUCT NUMBER 7703.10

Medical, health and well-Being



- › Approved Class III, MPG Class 1
- › Stable handrail for secure hold
- › Large platform 890 x 890 mm
- › Particularly stable
- › Low tread
- › Comfortable reading of weight value

ADIPOSIITY SCALE

PRODUCT NUMBER 7703.10



Characteristics

The Soehnle Professional Adiposity Scale is especially suited for dietary consultations, adiposity clinics, hospitals and rehabilitation institutes due to its large and low tread and the stable handrails. Special handholds at eye level provide orientation and safety to very adipose patients. In addition to EU Class III and the MPG certification, the high accuracy, the simple operation and the easy to read rotating indicator are further important factors for weighing operations with adipose, handicapped patients.

Version

- › Approved Class III, MPG Class 1
- › Platform size: 890 x 890 mm
- › Tread: 84 mm high
- › Indicator made of aluminium with plastic parts
- › 7-digit LCD-Indicator, back-lighted, digit height 20 mm
- › Indicator rotatable by 180° for more comfortable reading
- › BMI-, Hold- and dialysis function
- › Units in kg or lb

Technical data

- › Mains operation: power supply 230 V, 50/60 Hz, 2 m cable
- › Battery operation for approx. 40 hours with 4 Mignon-batteries 1,5 V (not included in delivery scope)
- › Auto-Off-function at battery operation (time adjustable)

Optional

- › Charger with rechargeable batteries for approx. 40 hours operation time (Order number 2557.04.001)
- › Serial Interface RS 232
- › Foldable height measuring rod (Order number 5003.01.001)

Order information

Order number	Weighing range	Division	Dimensions	Weight
	kg	g	mm	kg
7703.10.001 (without interface)	330	to 200 kg = 100 g, from 200 kg = 200 g	890 x 890 x 1484	35
7703.10.002 (with interface)	330	to 200 kg = 100 g, from 200 kg = 200 g	890 x 890 x 1484	35

Your Soehnle Professional dealer:

www.soehnle-professional.com