

# BLUETOOTH BED SCALE

PRODUCT NUMBER 7711

Medical, health and well-being



- › Approved to Class III, MPG Class 1
- › Wireless measured value transmission
- › Low power requirement of foot sensors
- › Range: 10m
- › EasyLift system - innovative hoist technology

# BLUETOOTH BED SCALE

PRODUCT NUMBER 7711



## Features

The handy bed scale with its functional design from Soehnle Professional is ideal for applications in hospitals, rehabilitation clinics and homes for the elderly. The four weighing elements interlinked by radio communication allow for very simple installation since it requires no inconvenient cables. This also enables much easier cleaning and ensures improved hygiene in hospital rooms. The evaluator includes optional interfaces to transfer weight values to further processing equipment.

### Technical data

- › Powered by battery or mains
- › Battery capacity 10 hrs on continuous duty, standby 125 hrs
- › 8 batteries included in scope of delivery (lead-gel rechargeable battery 6V, 1.2Ah)
- › Charger for 4 batteries (foot sensor) with quick charge and trickle charge
- › Wheel diameters 130–160mm (smaller diameters optional)
- › Wheel width max. 71mm
- › Interfaces to PC/network (optional)

### Design

- › Cordless measuring data transfer from foot sensor to digital indicator
- › Low power requirement of foot sensors
- › Lift facility integrated in foot sensor
- › Calm weight display through filtration
- › Range: 10m
- › Easy to read due to large digits
- › LCD, digit height: 20mm

### Options

- › Interface for transfer and further processing of measured values

## Ordering information

Order Number	Weighing Range	Division	Dimensions trolley	Weight
with trolley	kg	g	mm	kg
7711.01.001	500	up to 300kg = 100g, from 300kg = 200g	875 x 470 x 1 080	30.8
without trolley				
7711.02.001	500	up to 300kg = 100g, from 300kg = 200g		22.8

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)