

BABYSCALE

PRODUCT NUMBER 7725

Medical, health and well-being



- > Approved to Class III
- > MPG Class 1
- > Tare function
- > Hold function

BABYSCALE

PRODUCT NUMBER 7725

Features

The Soehnle Professional Babyscale 7725 is a two-range scale with high precision and reliability. The damping and hold function guarantee a stable weight indication. Another advantage of the babyscale is its high stability for safe weighing of our babies. Easy to clean plastic and ergonomic forms without chamfers and components allow for simple and hygienic cleaning.

- › Easy handling
- › High-contrast LCD display with 20 mm high figures
- › Large reading angle
- › Automatic zeroing at switch-on
- › Tare function key (laid on pad can be tared)
- › Hold function
- › Acoustic signals
- › Main and battery operation
- › Power plug with control LED (included in delivery scope)
- › Auto-Off at battery operation

Accessory

- › Baby height rod (Order No. 5005.01.001) to be mounted simply

Order information

Order Number	Weighing Range	Division	Dimensions	Weight
	kg	g	mm	kg
7725.01.001	15	up to 6 kg = 2 g, from 6 kg = 5 g	600 x 355 x 120	4.6

Your Soehnle Professional dealer:

www.soehnle-professional.com