

# WHEELCHAIR WITH INTEGRATED SCALE

PRODUCT NUMBER 7802.10

Medical, health & well-being



- › Weighing range: 120kg
- › Easy-running, foldable wheelchair
- › Simple and easy transport of the scale to the patient
- › Attachable indicator
- › Perfectly suited for caring institutions

# WHEELCHAIR WITH INTEGRATED SCALE

PRODUCT NUMBER 7802.10



## Characteristics

The Soehnle Professional Wheelchair is the perfect scale for weight control of patients with a mobility handicap. The personal weighing in a sitting position and the simple and easy transport of the scale is of particular advantage especially for elderly homes and caring institutions.

### Version

- › Scale with simple one-key operation
- › Zeroing and taring function
- › High contrast LCD display, digit height 20mm
- › Padded seat and backrest
- › 24 inch (600mm) wheels made of plastic with puncture-proof tyres
- › Foot rests hinged towards the inside and outside, adjustable, detachable
- › Armrest detachable
- › Handbrake
- › Foldable to 35cm
- › 4 load cells integrated
- › Indicator on pushbar detachable
- › Not intended for weight indication at practicing medical science

### Technical data

- › Mains and battery operation
- › Battery operation with 4 x 1.5V Mignon batteries for approx. 150 hours
- › Auto-off function in battery mode
- › Dimensions: 1,060 x 660 x 970mm
- › Seat height: 51cm, seat depth: 42cm, seat width: 46cm
- › Unit switchable to kg/lb

### Optional

- › Charger with rechargeable batteries for approx. 80 hours operation (Order number 2557.03.001)

## Order information

Order number	Version	Load capacity	Division	Weight
		kg	g	kg
7802.10.001	Battery operation, incl. mains adapter	120	50	22
7802.10.002	With rechargeable batteries, incl. charger	120	50	22

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)