

# BABY HEIGHT ROD

PRODUCT NUMBER 5001.10

Medical, health and well-being



# BABY HEIGHT ROD

PRODUCT NUMBER 5001.10



## Features

### Technical data

- > Easy to use, quick and precise
- > Slide made of white plastic
- > Rod made of metal
- > Easy to clean

## Ordering information

Order Numbers	Measuring Range	Divisions	Dimensions	Weight
	cm	cm	mm	g
5001.10.001	40 - 80	0.5	500 x 340	300

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)