

MEDISCALE ANALOGUE

PRODUCT NUMBER 6190

Medical, health and well-being



- › Approved to Class III
- › Large analogue indicator made of glass
- › Mechanical functioning
- › Weighing range 150kg

MEDISCALE ANALOGUE

PRODUCT NUMBER 6190



Features

The Soehnle Professional MediScale analogue is especially interesting for hospitals, surgeries and medical institutions that use approved medical products according to Class III. A large and easy-readable analogue scale, a low tread in stylish design, the easy use and the abdication of electronics recommend the scale's use in pharmacies, fitness studios and wellness locations.

Characteristics

- > Large analogue scale made of glass
- > Indication with a division of 1kg
- > Approved to Class III (not to be used according to MPG – Medical Product Law)
- > Units in kg or lbs
- > No batteries or power supply required
- > Large tread 370 x 325mm

Ordering information

Order Number	Weighing Range	Division	Dimensions	Weight
	kg	kg	mm	kg
6190.01.001	150	1	530 x 350 x 90	4.8

Your Soehnle Professional dealer:

www.soehnle-professional.com