

# PERSONAL SCALES

PRODUCT NUMBER 7700

Medical, health and well-being



- › Approved to EU Class III, MPG Class 1
- › BMI and dialysis functions
- › Weighing range 250kg
- › Divisions: from 50g
- › Mains and rechargeable battery power supplies

# PERSONAL SCALES

PRODUCT NUMBER 7700



## Features

Scales for ultra high demands. Besides approvals to EU Class III and MPG, the main factors for weighing in doctor's surgeries and hospitals are high precision, reliability, user-friendliness and good legibility.

### Technical data

- › Mains power supply: power adapter input 230V/50-60Hz, 2m power cord
- › Battery power supply for approx. 120 hrs with 4 x 1.5V Mignon cells
- › Auto-Off function for battery power supply (time adjustable)
- › Digital indicator dimensions: l x w x h 306 x 77 x 82.6mm
- › Platform size: l x w 310 x 275mm

### Design

- › Approved to EU Class III, MPG Class 1
- › BMI, Hold and dialysis functions
- › Weighing platform with easy to clean, skin-friendly covering
- › Digital indicator made of aluminium profile sections with plastic parts
- › With stand
- › 7-digit LCD, backlit, digit height 20mm
- › Two-range version
- › Units in kg or lb

### Options

- › Charge adapter with rechargeable batteries for approx. 120 hrs service time (Order No. 2557.04.001)
- › RS232 serial interface
- › Set of rollers for mobility 7700.01.001 (Order No. 2513.01.001)
- › Height rod (Order No. 5002.02.001)
- › Electronic body height measuring instrument (Order No. 5001.01.001)

## Ordering information

Order Numbers	Order Numbers	Weighing Range	Divisions	Dimensions	Weight
without interface	with interface	kg	g	mm	kg
7700.01.001	7700.01.002	250	up to 150kg = 50g over 150kg = 100g	310 x 400 x 865	12

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)