

PERSONAL SCALES

PRODUCT NUMBER 7701

Medical, health and well-being



- › Approved to EU Class III, MPG Class 1
- › BMI and dialysis functions
- › Weighing range 250kg
- › Divisions: from 50g
- › Mains and rechargeable battery power supplies

PERSONAL SCALES

PRODUCT NUMBER 7701



Features

Scales for ultra high demands – Besides approvals to EU Class III and MPG, the main factors for weighing in doctor's surgeries and hospitals are high precision, reliability, user-friendliness and good legibility.

Technical data

- > Mains power supply: power adapter input 230V/50-60Hz, 2m power cord
- > Battery power supply for approx. 160 hours with 4 x 1.5V Mignon cells (not contained in scope of delivery)
- > Auto-Off function for battery power supply (time adjustable)
- > Digital indicator dimensions: l x w x h 306 x 77 x 82.6mm
- > Platform size: l x w 310 x 275mm
- > The measuring cable between weighing platform and digital indicator is approx. 2m long.

Design

- > Approved to EU Class III, MPG Class 1
- > BMI and dialysis functions

- > Weighing platform with easy to clean, skin-friendly covering
- > Digital indicator made of aluminium profile sections with plastic parts and bench or wall bracket
- > 7-digit LCD, backlit, digit height 20mm
- > Two-range version
- > Units in kg or lb

Options

- > Charge adapter with rechargeable batteries for approx. 120 hrs service time (Order No. 2557.04.001)
- > RS 232 serial interface
- > Electronic body height measuring instrument (Order No. 5001.01.001)
- > With carry case

Ordering information

Order Numbers	Weighing Range	Divisions	Dimensions	Weight
without interface	kg	g	mm	kg
without carry case:				
7701.02.001	250	up to 150kg = 50g, from 150kg = 100g	400 x 310	8
with carry case:				
7701.01.001	250	up to 150kg = 50g, from 150kg = 100g	605 x 395 x 135	10
with interface				
without carry case:				
7701.02.002	250	up to 150kg = 50g, from 150kg = 100g	400 x 310	8
with carry case:				
7701.01.002	250	up to 150kg = 50g, from 150kg = 100g	605 x 395 x 135	10

Your Soehnle Professional dealer:

www.soehnle-professional.com