

# ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7702

Medical, health and well-being



- › Electronic chair weigher
- › Approved to Class III, MPG Class 1
- › Hinged footrests and armrests
- › No wear parts and maintenance-free
- › Weighing range 200kg
- › BMI and dialysis functions

# ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7702



## Features

Electronic chair weighers to meet very high requirements – besides approval to EU Class III and to MPG, the main features include high precision and reliability, user-friendliness and easy readability. These features are ideally suited for hospitals, old people's homes and nursing homes for weighing patients who have difficulty walking or standing.

### Design

- > Approved to EU Class III, MPG Class 1
- > Seat made of easy to clean plastic, ergonomically shaped
- > Hinged footrests and armrests
- > 4 wheels, of which 2 lockable steering wheels
- > Digital indicator: housing made of aluminium sections with plastic parts
- > 7-digit LCD, backlit, digit height 20mm
- > Two-range version
- > Units in kg or lb

### Technical data

- > Mains power: power adapter input 230VAC, 50Hz, output 9VDC, 300mA
- > Battery power supply for approx. 160 hrs by 4 x 1.5V Mignon cells
- > Auto-Off function (time adjustable)
- > Dimensions: lxxh 580 x 1100 x 930mm
- > Unit switchover kg/lb

### Options

- > Charge adapter with rechargeable batteries for approx. 120 hrs service time (Order No. 2557.04.001)
- > Pouch for power adapter for attachment to wheelchair scale (Order No. 2507.03.001)
- > RS 232 serial interface for connection to PC/printer

## Ordering information

Order Numbers	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7702.01.001 (without interface)	200	up to 150kg = 50g over 150kg = 100g	580 x 1100 x 930	27
7702.01.002 (with interface)	200	up to 150kg = 50g over 150kg = 100g	580 x 1100 x 930	27

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)