


WHEELCHAIR WEIGHER

PRODUCT NUMBER 7705.10

Medical, health and well-being



- › Approved to Class III, MPG Class 1
- › Transportable by means of 2 rollers and handle
- › Large weighing platform
- › Mains power supply and rechargeable battery

WHEELCHAIR WEIGHER

PRODUCT NUMBER 7705.10



Features

The Electronic Wheelchair Weigher 7705 is the ideal scale to check the weight of patients in wheelchairs. Besides approval to Class III and MPG approval, precision and reliability, user-friendliness and clear legibility are the main factors for weight control in hospitals, old people's homes and nursing homes.

Technical data

- > Mains power: power adapter 230V/50-60Hz, 2m cable
- > Battery power supply for approx. 40 hrs by 4 x 1.5V Mignon cells (not contained in the scope of delivery)
- > Auto-Off function in battery power mode (time adjustable)
- > Dimensions: l x h x d 950 x 800 x 55mm

Design

- > Approved to Class III, MPG Class 1
- > Footplate: 55mm high
- > Transportable by means of 2 rollers and handle
- > Indicator made of aluminium profile section with plastic parts
- > 7-digit LCD, backlit, digit height 20mm
- > Taring function: wheelchairs and clothing are tarable
- > BMI, Hold and dialysis functions
- > Units in kg or lb

Options

- > Charge adapter with rechargeable batteries for approx. 40 hrs service time (Order No. 2557.04.001)
- > RS 232 serial interface
- > Electronic body height measuring instrument (Order No. 5001.01.001)

Ordering information

Order Numbers	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7705.10.001 (without interface)	300	100	950 x 800 x 55	28
7705.10.002 (with interface)	300	100	950 x 800 x 55	28

Your Soehnle Professional dealer:

www.soehnle-professional.com