

# LIFTER SCALES

PRODUCT NUMBER 7707

Medical, health and well-being



- › Approved to EU Class III, MPG Class 1
- › Integratable in many lifter systems
- › Designed for retrofitting
- › Weighing range 250kg
- › Mains and rechargeable battery power supplies

# LIFTER SCALES

PRODUCT NUMBER 7707



## Features

Electronic lifter scales to meet ultra high demands – Besides approvals to EU Class III and MPG, the main factors for weighing bedridden patients in hospitals, residential and nursing homes are high precision, reliability, user-friendliness and easy legibility.

### Technical data

- › Mains power supply: power adapter 230V/50-60Hz, 2m power cord
- › Battery power supply with 4 x 1.5V Mignon cells (not contained in scope of delivery)
- › Auto-Off function for battery power supply (time adjustable)
- › Backlit LCD with 20mm high digits
- › Weighing element dimensions: lxwxh 200 x 70 x 25mm

### Design

- › Approved to EU Class III, MPG Class 1
- › BMI, Hold and dialysis functions
- › Easy to use and sturdy design
- › Digital indicator made of aluminium profile sections with plastic parts
- › 7-digit LCD, backlit, digit height 20mm
- › Units in kg or lb

### Options

- › Charge adapter with rechargeable batteries for approx. 120 hrs service time (Order No. 2557.04.001)
- › RS 232 serial interface

## Ordering information

Order Numbers	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7707.01.001 (without interface)	250	100	200 X 70 X 25	2
7707.01.002 (with interface)	250	100	200 X 70 X 25	2

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)