

# BED SCALES

PRODUCT NUMBER 7710

Medical, health and well-being



- › Approved to EU Class III, MPG Class 1
- › EasyLift system - innovative hoist technology
- › Alphanumeric input keypad
- › Weighing range 500kg
- › Divisions: 100g/200g
- › No wear parts and maintenance-free

# BED SCALES

PRODUCT NUMBER 7710



## Features

The scales for weight control in intensive care units and dialysis centres, consisting of 4 weighing units, Digital Indicator 3020 and a practical trolley.

- › New simplified handling with EasyLift system – innovative hoist technology
- › No more rolling beds onto weighing facilities. This is especially good for patients under intensive pain
- › The weighing units are simply positioned under the bed wheels from the side. The bed is raised by foot pedal without the use of excessive force.
- › The weighing units are individually plugged in

### Technical data

- › Mains power supply: power adapter 100-240V/50-60 Hz, 400mA, 2m power cord
- › 4 weighing units, digital indicator made of stainless steel with 2m power cord
- › Units in kg or lb

### Design

- › BMI, Hold and dialysis functions
- › Alphanumeric keypad to enter organisation data
- › RS 232 interface to transfer data to IT network or printer
- › Plus/minus weighing
- › Approved to MPG Class 1 and EU Class III
- › Dot matrix display with 240 x 64 pixels
- › Membrane keyboard with large keys and pressure point

### Options

- › Tape printer (Order No. 2795.11.001)
- › Data cable (Order No. 2550.03.002)
- › Trolley

## Ordering information

Order Numbers	Weighing Range	Dimensions	Weight	Divisions
with trolley	kg	mm	kg	g
7710.01.001	500	860 x 480 x 1050	33	up to 300kg = 100g, from 300kg = 200g
without trolley				
7710.02.001	500		25	up to 300kg = 100g, from 300kg = 200g

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)