

# BABY SCALE

PRODUCT NUMBER 7726.03

Medical, health and well-being



- > Approval Class III
- > MPG Class 1
- > Backlit LCD with 20mm high digits
- > Hold function

# BABY SCALE

PRODUCT NUMBER 7726.03



## Features

Two-range scale with high precision and reliability, stable weight display by damping and hold function, high stability baby weighing tray made of easy to clean plastic.

- › Easy to operate
- › High contrast, backlit LCD with 20mm high digits
- › Automatic zeroing on power-on
- › Tare button (to tare blankets, etc.)
- › Print and data transfer button
- › RS 232 data interface (V24)
- › Mains power supply or rechargeable battery
- › With weight lock function
- › Unit indicator switchover kg/lb

### Accessories

- › Baby height rod (Order No. 5001.10.001)

## Ordering information

Order Numbers	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7726.03.001	15	up to 6kg = 2g from 6kg = 5g	530 x 320 x 180	6.6

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)