


VETERINARY SCALES

PRODUCT NUMBER 7758.10

Medical, health and well-being



- › Transportable via a handle
- › Light-weight construction
- › Large platform
- › Extra flat, only 35 mm high
- › Platform made of stainless steel
- › Mains and battery operation

VETERINARY SCALES

PRODUCT NUMBER 7758.10



Characteristics

The electronic Veterinary Scales 7758 are the ideal scales for weight control of animals in sectors, which are not subject to the Medical Device regulation. Uncomplicated and mobile handling, robust technology, large indicator and ergonomic, modern design characterise Veterinary Scales of Soehnle Professional.

Technical data

- › Mains power supply: power adapter 230V/50-60 Hz, 10mA, 2.2VA, 2m power cable
- › Battery power supply with 4 x 1.5V Mignon cells (not included in scope of delivery)
- › Auto-off function in battery mode
- › Platform size: L 900 x W 600 mm

Functions

- › Switch on/off, hold function

Design

- › Easy to use and sturdy design
- › Ultra flat design 35mm high
- › With black rubber mat
- › Transportable via a handle
- › Digital indicator: housing made of impact resistant plastic
- › High-contrast LCD with 20mm high digits
- › Simple one-key operation
- › Unit switchover kg/lb

Options

- › Charge adapter with rechargeable batteries for approx. 150 hrs service time (Order No. 2557.03.001)

Ordering information

Order numbers	Weighing range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7758.10.001 (with batteries) incl. mains adapter	200	100	900 x 600 x 35	12
7758.10.002 (with rechargeable batteries) incl. charger	200	100	900 x 600 x 35	12

Your Soehnle Professional dealer:

www.soehnle-professional.com