

VETERINARY SCALES

PRODUCT NUMBER 7758

Medical, health and well-being



- › Transportable by means of 2 rollers and handle
- › Large platform 900 x 600mm
- › Ultra flat design 55mm high
- › Platform made of chrome nickel steel
- › Mains and rechargeable battery power supplies

VETERINARY SCALES

PRODUCT NUMBER 7758



Features

The electronic Veterinary Scales 7758 are the ideal weight control scales for animals in sectors which are not subject to the Medical Device Ordinance. Soehnle Professional veterinary scales are easy to use and mobile, have a sturdy design, large display and an ergonomic, modern design.

Technical data

- › Mains power supply: power adapter 230V, 50–60 Hz, 10mA, 2.2VA, 2m power cord
- › Battery power supply with 4 x 1.5V Mignon cells (not contained in scope of delivery)
- › Auto-Off function in battery power mode
- › Platform size: lxxwx 900 x 600 x 55mm

Design

- › Easy to use and sturdy design
- › Ultra flat design 55mm high
- › With black rubber mat
- › Transportable by means of 2 rollers and handle
- › Digital indicator: housing made of impact resistant plastic
- › High-contrast LCD with 20mm high digits
- › Simple one-key operation
- › Unit switchover kg/lb

Options

- › Charge adapter with rechargeable batteries for approx. 150 hrs service time (Order No. 2557.03.001)

Ordering information

Order Number	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7758.01.001 (powered by battery cells) incl. power adapter	200	100	900 x 600 x 55	22
7758.01.002 (powered by rechargeable batteries) incl. 2557.03.001	200	100	900 x 600 x 55	22

Your Soehnle Professional dealer:

www.soehnle-professional.com