

# PERSONAL SCALES

PRODUCT NUMBER 7800

Medical, health and well-being



- › Mains and rechargeable battery power supplies
- › No wear parts and maintenance-free
- › Auto-Off function in battery power mode

# PERSONAL SCALES

PRODUCT NUMBER 7800



## Features

The electronic Personal Scales 7800 are the ideal weight control scales in sectors which are not subject to the Medical Device Ordinance. Soehnle Professional personal scales are easy to use, have a sturdy design, large display and an ergonomic, modern design.

### Technical data

- > Mains power supply: power adapter 230V/50/60Hz, 10mA, 2.2VA, 2m power cord
- > Battery power supply for approx. 350 hours with 4 x 1.5V Mignon cells (not contained in scope of delivery)
- > Auto-Off function in battery power mode
- > Platform size: l x w 310 x 275mm

### Design

- > QUATTRO-TRONIC: platform with 4 load cells
- > Weighing platform with easy to clean, skin-friendly covering
- > Ultra flat platform

- > Digital indicator: housing made of impact resistant plastic
- > No wear parts and maintenance-free
- > Simple one-key operation
- > With stand
- > Unit switchover kg/lb

### Options

- > Charge adapter with rechargeable batteries for approx. 200 hrs continuous service time (Order No. 2557.03.001)
- > Height rod (Order No. 5002.02.001)
- > Electronic body height measuring instrument (Order No. 5001.01.001)

## Ordering information

Order Numbers	Order Numbers	Weighing Range	Divisions	Dimensions	Weight
powered by batteries	with rechargeable battery	kg	g	mm	kg
7800.01.001	7800.01.002	200	100	450 x 350 x 860	10
incl. power adapter	incl. 2557.03.001				

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)