

ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7802

Medical, health and well-being



- › Weighing range 200kg
- › Hinged footrests and armrests
- › No wear parts and maintenance-free

ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7802



Features

The mobile Electronic Chair Weigher 7802 is the ideal scale for checking the weight of people who have difficulty walking or standing in areas not subject to the Medical Product Law. The Soehnle Professional electronic chair weigher is easy to use and has a sturdy design, large indicator, comfortable seat and ergonomic timeless design.

Design

- > Seat made of easy to clean plastic, ergonomically shaped
- > Hinged footrests
- > Hinged armrests
- > 4 wheels, of which 2 lockable steering wheels
- > Digital indicator: housing made of impact resistant plastic
- > No wear parts and maintenance-free
- > Simple one-key operation

Technical data

- > Mains power data: 230V, 50/60Hz, 10mA, 2.2VA
- > Battery power supply for approx. 150 hrs with 4 x 1.5V Mignon cells
- > Auto-Off function in battery power mode
- > Dimensions: l x w x h 580 x 1100 x 930mm
- > Unit switchover kg/lb

Options

- > Charge adapter with rechargeable batteries for approx. 80 hrs continuous service time (Order No. 2557.03.001)

Ordering information

Order Numbers	Order Numbers	Weighing Range	Division	Dimensions	Weight
rechargeable battery	rechargeable battery	kg	g	mm	kg
mains power supply	incl. charge adapter				
7802.01.001	7802.01.002	200	100	580 x 1100 x 930	27
incl. power supply	incl. 2557.03.001				

Your Soehnle Professional dealer:

www.soehnle-professional.com