




7803.50

Multi-function scale

for adipose patients and patients in wheelchairs




- S** Can be freely combined and retrofitted
- S** Stable handrails for secure hold

*Medical, health and
well-being*

Multi-function scale

for adipose patients and patients in wheelchairs



Weighing platform with two rollers and a handle, incl. indicator

- › Transportable by means of two rollers and handle
- › Handrail, drive-on ramps, side protection and adiposity handrail freely combinable
- › Weighing platform and handrails made of aluminium
- › Large platform (800 x 800 mm with handrails or side protection)
- › Secure hold with sturdy handrails
- › Indicator can be rotated by 180°
- › Easy-to-operate indicator
- › Backlit LCD display, digit height 20 mm
- › Wall mounting/bench bracket for indicator included in scope of delivery

Order information

| | |
|----------------|-------------------|
| Order number | 7803.50.001 |
| Weighing range | 300 kg |
| Division | 100 g |
| Dimensions | 840 x 880 x 55 mm |

Combination options/accessories

The following handrails, drive-on ramps and side guides can be fitted to approved scale

Support handrails

The following order number comprises 1 handrail. Up to 3 handrails can be fitted. (By ordering two support rails per scale, 1 connecting element for the railing is included in scope of delivery. By ordering 3 support rails two connecting elements are included in scope of delivery)


Order number 2530.03.001
Dimensions 750 x 850 x 30 mm

Adiposity handrails (set)

The following order number comprises 2 handrails and 1 support handrail 2530.03.001.


(2 connecting elements for the railing is included in scope of delivery)

Order number 2530.04.001
Dimensions 1,360 x 850 x 30 mm




Drive-on ramp designed to accommodate wheelchairs. The following order number comprises 1 ramp. 2 ramps can be fitted.

Order number 2505.08.001
Dimensions 880 x 75 x 48 mm



Side guide (set) for use as a wheelchair weigher. The following order number comprises 2 side profiles.


Order number 2502.04.001
Dimensions 30 x 30 x 820 mm



Handrail connection

Connecting element for a better stability of the handrails. The following order number comprises 1 handrail connection, with wrench.

Order number 2510.07.001



Handrail bracket for indicator

To mount the indicator on the support handrails of the scale.

Order number 2510.08.002



Charger with rechargeable batteries

Approx. 200 hours of operation.

Order number 2557.03.001

www.soehnle-professional.com

Your Soehnle Professional dealer: