


ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7805

Medical, health and well-being



ELECTRONIC CHAIR WEIGHER

PRODUCT NUMBER 7805



Features

The Electronic Wheelchair Weigher 7805 is the ideal scale to check the weight of people in wheelchairs who have difficulty walking or standing. Possible applications in pharmacies, rehabilitation clinics or old people's homes. The Soehnle Professional electronic wheelchair weigher features ease of use, robust electronics, ultra flat platform with integrated ramps, a large indicator and ergonomic timeless design.

Design

- › Transportable by means of 2 rolls and handle
- › Extremely lightweight platform weighing only 14kg, aluminium frame, weighing platform made of waterproof bonded laminated wood
- › Non-slip surface
- › Digital indicator: housing made of impact resistant plastic, with 1.5m cable
- › Wall bracket for digital indicator included in scope of delivery
- › No wear parts and maintenance-free
- › Simple one-key operation

Technical data

- › Mains powered: 230V, 50/60Hz
- › Battery power supply for approx. 150 hrs with 4 x 1.5V Mignon cells
- › Auto-Off function in battery power mode
- › Unit switchover kg/lb

Options

- › Charge adapter with rechargeable batteries for approx. 80 hrs continuous service time (Order No. 2557.03.001)

Ordering information

Order Numbers	Order Numbers	Weighing Range	Division	Dimensions	Weight
rechargeable battery	rechargeable battery	kg	g	mm	kg
mains power supply	incl. charge adapter				
7805.01.001	7805.01.002	200	100	800 x 800 x 40	14
incl. power supply	incl. 2557.03.001				

Your Soehnle Professional dealer:

www.soehnle-professional.com