

# ADIPOSIITY SCALES

PRODUCT NUMBER 7806

Medical, health and well-being



- › Weighing range 300kg
- › Large platform 600 x 600mm
- › Ultra flat design 40mm high
- › With handle
- › Mains and rechargeable battery power supplies

# ADIPOSIITY SCALES

PRODUCT NUMBER 7806



## Features

The electronic Adiposity Scales 7806 are the ideal weight control scales for adiposity patients in sectors which are not subject to the Medical Device Ordinance. Soehnle Professional adiposity scales are easy to use, have a sturdy design, large display and an ergonomic, modern design.

### Technical data

- › Mains power supply: power adapter 230V/50-60 Hz, 10mA, 2.2VA, 2m power cord
- › Battery power supply for approx. 150 hours with 4 x 1.5V Mignon cells (not contained in scope of delivery)
- › Auto-Off function in battery power mode
- › Dimensions: lxxwxh 600 x 600 x 40mm

### Design

- › Simple one-key operation
- › Easy to use and sturdy design
- › Ultra flat design 40mm high.
- › With handle
- › Digital indicator: housing made of impact resistant plastic
- › High-contrast LCD with 20mm high digits
- › Unit switchover kg/lb

### Options

- › Charge adapter with rechargeable batteries for approx. 150 hrs service time  
(Order No. 2557.03.001)
- › Electronic body height measuring instrument  
(Order No. 5001.01.001)

## Ordering information

Order Numbers	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
7806.01.001 (battery-powered) incl. power adapter	300	100	600 x 600 x 40	18
7806.01.002 (powered by rechargeable batteries) incl. 2557.03.001	300	100	600 x 600 x 40	18

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)