

MEDISCALE DIGITAL

PRODUCT NUMBER 7840

Medical, health and well-being



- › Exclusive personal scale with digital and analogue indicator
- › Large tread with antislip surface
- › Weighing range 180kg
- › Division 100g

MEDISCALE DIGITAL

PRODUCT NUMBER 7840

Characteristics

The fully electronic MediScale is especially suitable for usage in pharmacies, fitness and wellness institutions as well as for nutrition counselling. With an extra flat and big weighing platform, the weighing is also easy for insecure people. The so far unknown combination of the familiar analogue indication with the exact readability of the digital indication makes the scale to a special highlight.

Version

- > Personal scale with digital and analogue glass indicator
- > Immediate start function
- > Large tread 405 x 325mm
- > Antislip surface
- > LCD-Display with 25mm high digits
- > Powered operated analogue indication

Technical data

- > Battery operation with 4 x 1.5V AA alkaline batteries
(not included in delivery scope)
- > Auto-Off function

Order information

Order Number	Weighing Range	Division	Dimensions	Weight
	kg	g	mm	kg
7840.01.001	180	100	550 X 325 X 120	5.4

Your Soehnle Professional dealer:

www.soehnle-professional.com