

# ANALYTICAL SCALE STANDARD

PRODUCT NUMBER 9401

Laboratory scales for medical and industrial applications



- > Weighing range: 160 - 310 g
- > Division 0.1 mg
- > Large graphical and back-lighted display
- > Automatic internal calibration

# ANALYTICAL SCALE STANDARD

PRODUCT NUMBER 9401

## Features



The series 9401 is the middle class of Analytical Scales. It has a back-lighted graphical display as well as an automatic internal calibration, which is responsive to temperature changes. Or also works time-controlled. Numerous functions make different applications possible.

### Functions

- › Piece counting
- › Percentage weighing
- › Checking
- › Dosing
- › Recipes
- › Animal weighing
- › Density determination
- › Statistics
- › Air buoyancy correction
- › GLP procedures
- › Underfloor weighing
- › Programmable printouts

### Technical data

- › Operating temperature: + 10° till + 45° C
- › Power supply: Mains adapter 230 V 50 Hz AC, 9 V DC
- › Automatic internal calibration
- › Menu navigation: 6 languages selectable
- › Stabilisation time: 3 - 5 second
- › Interface RS 232
- › PS2 interface for keyboard
- › Real time clock
- › Power supply included in delivery scope
- › Interface for an additional display possible

### Optional

- › Density determination set
- › Platform for underfloor weighing
- › Thermal printer with clock module
- › Thermal printer with statistics function
- › Tape printer
- › Anti-vibration desk made of steel
- › PC software for transmitting and saving measuring data

## Ordering information

Order number	Weighing range	Division		Dimensions	Platform size
	g	mg (d)	mg (e)	mm	Diameter mm
9401.01.005	160	0.1	1	335 X 210 X 335	85
9401.01.010	220	0.1	1	335 X 210 X 335	85
9401.01.020	310	0.1	1	335 X 210 X 335	85

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)