

# ANALYTICAL SCALE PRO

PRODUCT NUMBER 9402

Laboratory scales for medical and industrial applications



- > Weighing range: 60 – 310 g
- > Division 0.01 g
- > Large graphical and back-lighted display, removable
- > Automatic internal calibration
- > Approved class I

# ANALYTICAL SCALE PRO

PRODUCT NUMBER 9402

## Features



The series 9402 is the luxury class of Analytical Scales. It has a back-lighted graphical display as well as an automatic internal calibration, which is responsive to temperature changes or also works time-controlled. and is perfect for fast definition of measure in laboratory- and industrial sector. Beyond the high resolution and the possibility to detach the display characterizes the scale.

### Functions

- › Piece counting
- › Percentage weighing
- › Checking
- › Dosing
- › Recipes
- › Animal weighing
- › Density determination
- › Pipette calibration
- › Statistics
- › Air buoyancy correction
- › GLP procedures
- › Underfloor weighing
- › Programmable printouts

### Technical data

- › Principle of measurement: Electromagnetic power compensation
- › Operating temperature: + 10° till + 50° C
- › Power supply: Mains adapter 230 V 50 Hz AC, 9 V DC
- › Menu navigation: 6 languages selectable
- › Automatic internal calibration
- › Stabilisation time: 3 – 5 second
- › Interface RS 232
- › PS2 interface for keyboard
- › Real time clock
- › Power supply included in delivery scope
- › Interface for an additional display possible

### Optional

- › Density determination set
- › Pipettes calibration kit
- › Platform for underfloor weighing
- › Thermal printer with clock module
- › Thermal printer with statistics function
- › Tape printer
- › Anti-vibration desk made of steel
- › PC software for transmitting and saving measuring data

## Ordering information

Order number	Weighing range		Division		Dimensions	Platform size
	g	mg (d)	mg (e)	mm	Diameter mm	
9402.01.010	100	0.01	1	484 x 202 x 294	85	
9402.01.020	210	0.01	1	484 x 202 x 294	85	
9402.01.030	60 / 220	0.01 / 0.1	1	484 x 202 x 294	85	
9402.01.040	220	0.1	1	484 x 202 x 294	100	
9402.01.050	310	0.1	1	484 x 202 x 294	100	

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)