


# PRECISION SCALE STANDARD

PRODUCT NUMBER 9431

Laboratory scales for medical and industrial applications



- > Weighing range: 200 g – 6,000 g
- > Division 0.001 g
- > Back-lighted LCD display
- > Automatic internal calibration

# PRECISION SCALE STANDARD

PRODUCT NUMBER 9431

## Features



The series 9431 is the middle class of Precision Scales and is perfect for fast definition of measure in laboratory- and industrial sector. It has a weighing pan of stainless steel, a back-lighted graphical display and marked by a high resolution.

### Functions

- › Piece counting
- › Percentage weighing
- › Totalising
- › Checking
- › Density determination\*
- › GLP procedures
- › Underfloor weighing

### Technical data

- › Operating temperature: + 15° till + 30° C
- › Power supply: Mains adapter 230 V AC, 9 V DC
- › Automatic internal calibration
- › Stabilisation time: 2 seconds
- › Interface RS 232
- › Power supply included in delivery scope
- › Interface for an additional display possible

### Optional

- › Density determination set
- › Platform for underfloor weighing
- › Thermal printer with clock module
- › Thermal printer with statistics function
- › Tape printer
- › Wind deflector
- › Anti-vibration desk made of steel
- › PC software for transmitting and saving measuring data

## Ordering information

| Order number | Weighing range | Division     |            | Dimensions     | Platform size |
|--------------|----------------|--------------|------------|----------------|---------------|
|              | g              | g (d)        | g (e)      | mm             | mm            |
| 9431.01.010* | 600            | 0.001        | 0.01       | 335 x 210 x 80 | 128 x 128     |
| 9431.01.020* | 200 / 2,000    | 0.001 / 0.01 | 0.01 / 0.1 | 335 x 210 x 80 | 128 x 128     |
| 9431.01.030  | 2,100          | 0.01         | 0.1        | 335 x 210 x 80 | 195 x 195     |
| 9431.01.040  | 4,500          | 0.01         | 0.1        | 335 x 210 x 80 | 195 x 195     |
| 9431.01.050  | 6,000          | 0.01         | 0.1        | 335 x 210 x 80 | 195 x 195     |

Your Soehnle Professional dealer:

[www.soehnle-professional.com](http://www.soehnle-professional.com)