SOEHNLE PROFESSIONAL



Baby scale "Neo" conformity assessed

6752

- Model 6752: Conformity assessed according to class III, MPG class 1 m
- Very fine division 5 g
- Easy-to-read LC display with 20 mm high digits and backlight
- Integrated length scale in the weighing tray
- Battery and mains operation







Features

- Model 6752: Conformity assessed according to Class III, MPG Class 1 m
- Particularly small division 5 g
- Easy-to-read LC display with 20 mm high digits and backlight
- Tare function: documents placed on the scale can be tared
- Tare hand entry possible via the tare button
- Hold function for easy reading of the weight value after weighing
- Integrated RS232 and USB interfaces for connecting a PC
- Integrated length scale in the weighing pan (units in 2.5 mm increments)
- Weighing pan material: Acrylic
- Units in kg
- Energy-saving auto-off function
- Battery and mains operation 100 - 240 V, 50-60 Hz (not included in scope of delivery)
- Battery status indicator

Products

Order number Variation	**Approvable**	**Division** g	**Maximum load** kg
6752.01.002 Baby scale "Neo" conformity assessed	Yes	5	15

Accessories

Order number Variation		Description
2550.03.022		USB data cable, cable length 3 m
5005.11.003		Baby height rod